Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

MUAC, Child 11.5 Red

Product Code :

SLE/AE/07



Science Lab
Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

MUAC, Child 11.5 Red (Product Code: SLE/AE/07) is a colour-banded tape used to measure a child's mid-upper arm circumference during nutrition screening by community health workers and clinic staff. The reading falls in a red zone up to 11.5 cm, a yellow zone from 11.5 to 12.5 cm or a green zone above, and the tape is also graduated in 1 mm steps up to 26.5 cm. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a lab instruments supplier to health programmes.

Specifications

Parameter	Detail
Type	Child MUAC measuring tape, colour coded
Range	Up to 26.5 cm
Graduation	1 mm, with a thicker line at 21.0 cm
Accuracy	±1 mm of the maximum measurement (26.5 cm)
Red zone	0 to 11.5 cm (Pantone 1795 C)
Yellow zone	11.5 to 12.5 cm (Pantone 107 C)
Green zone	12.5 to 26.5 cm (Pantone 369 C)
Operating temperature	0°C to +40°C
Material	Non-tear, stretch-resistant plasticised synthetic paper; minimum thickness 0.3 mm
Print	Permanent, solvent resistant, easily readable in low light
Product code	SLE/AE/07
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Child MUAC Tape Used For?

The child's arm hangs loosely while the tape is placed at the midpoint between shoulder and elbow. It is passed around the arm and pulled just snug, not so loose that it gapes and not so tight that it pinches the skin, and both the colour under the reading point and the millimetre value are noted. How each colour is acted on, and for which age group, follows the national or programme protocol the team works under.

Because health workers and volunteers can learn the colour method quickly, programmes preparing laboratory equipment tenders in India or abroad usually order these tapes in large numbers for door-to-door and camp screening. Keep a spare tape in each kit, as a creased or cut tape should be withdrawn from use.

What to Specify When Ordering

-
- her
-