Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Microscope Case

Product Code :

SLE/M/62



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Microscope Case (Product Code: SLE/M/62) is a hard carry case with a padded interior, used to store a microscope and move it safely between laboratories, classrooms, exam centres and outreach camps. Wide straps hold the instrument steady, inside pockets keep eyepieces, objective cases and manuals with it, and rubber feet allow the cases to be stacked. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, among its lab equipment accessories for microscopy.

Specifications

Parameter	Detail
Case type	Hard carry case
Interior	Padded, with wide straps to hold the microscope securely
Storage pockets	Additional inside pockets for eyepieces, objective cases, manuals and other items
Base	Rubber feet
Stacking	Stackable
Inside dimensions	16 in H x 7 in W x 10 in D
Product code	SLE/M/62
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Microscope Case Used For?

Microscopes are most often damaged not during use but while being carried or stored: a knocked stage, a scratched objective or a lost eyepiece can take an instrument out of service for a whole term. A hard, padded case keeps dust off between practicals, protects the optics when microscopes travel to exam halls, science fairs or rural outreach camps, and keeps small accessories with the instrument they belong to. Stackable cases with rubber feet also make store-room shelving tidier and steadier.

Before packing, remove any slide, rack the stage down, turn the lowest-power objective into position and switch off and cool any lamp. Secure the microscope with the straps, put spare eyepieces and objective cases in the pockets, and store the case upright. Measure your microscope first, including knobs, mirror or lamp housing, to be sure it sits within the inside dimensions.

What to Specify When Ordering

- Height, width and depth of the microscope to be carried, so fit can be checked against 16 x 7 x 10 in
- Number of cases required, one per microscope
- Any slide cases or cleaning items to include on the same quotation

Sup

This d the