Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Microscope Slide Warming Table

Product Code :

SLE/M/80



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Microscope Slide Warming Table (Product Code: SLE/M/80) is a thermostatically controlled warming plate for glass slides, used to dry smears, settle and adhere tissue sections and warm slides before mounting in histology, pathology-training and biology laboratories. Its 25 x 8 in black anodised aluminium surface holds fifty-eight 1 x 3 in slides at temperatures from ambient to 70°C, read on an LED display. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, as a biology lab equipment supplier in India.

Specifications

Parameter	Detail
Working surface	25 x 8 in, black anodised aluminium, gives contrast for slides
Capacity	Fifty-eight 1 x 3 in slides
Temperature control	Thermostatic
Temperature range	Ambient to 70°C
Display	LED temperature display
Indicator	Pilot light
Power	120 V AC, with power cord
Product code	SLE/M/80
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Slide Warming Table Used For?

In section cutting, paraffin ribbons from the microtome are floated out, picked up on slides and then need gentle, even heat to flatten, dry and stick firmly before staining. A warming table does that for a whole batch at once instead of one slide at a time over a flame. The same surface speeds the drying of blood and bacterial smears, warms slides before mounting and helps mounting media set. The black anodised finish makes clear glass and thin, pale sections easy to see and position.

Switch on and check the pilot light, set the thermostat to the temperature your protocol calls for within the ambient to 70°C range, and wait for the LED display to settle before loading slides. Lay slides flat in rows so none overhang the edge. The surface stays hot after use, so let it cool before wiping off wax or stain, and keep flammable clearing solvents away from the unit. Because it is listed for 120 V AC, confirm the supply version for your site when enquiring.

What to Specify When Ordering

- Mains supply at your site, since the listed model runs on 120 V AC
- Number of warming tables, allowing fifty-eight 1 x 3 in slides per table
- Slides, drying racks or a microtome to include on the same quotation

Sup

Jain S Lab
name