



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Mixer Artery Forceps

**Product Code :**

SLE/SI/30



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Mixer Artery Forceps** (Product Code: SLE/SI/30) is a springy threading forceps used to dissect around a vein or artery and pass a ligature thread behind it during surgery. Its 23 cm length and jaws that grip thread well make it a practical choice when a vessel has to be tied in a deep field. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, to hospital buyers at home and through lab equipment export from India to international markets.

## Specifications

| Parameter    | Detail                                                         |
|--------------|----------------------------------------------------------------|
| Type         | Threading forceps, springy                                     |
| Jaws         | Grip the thread well                                           |
| Ratchet      | Flexible, lockable                                             |
| Length       | 23 cm                                                          |
| Material     | Martensitic steel (quenched, magnetic)                         |
| Composition  | 0.20% carbon; 13% chromium                                     |
| Intended use | Dissecting; passing a thread around a vein or artery           |
| Product code | SLE/SI/30                                                      |
| Supplied by  | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India |

## What Is a Mixer Forceps Used For?

Before a vessel can be tied it has to be freed from the surrounding tissue on all sides. The surgeon works the closed tips of the Mixer forceps behind the vessel, opens them gently to separate the tissue, then grasps a ligature thread and draws it back around the vein or artery so that it can be tied.

The 23 cm length lets this be done deep in the operative field, and the lockable ratchet keeps the thread held while it is passed. Surgical skills labs that buy laboratory instruments in India can pair it with a long haemostat for ligation practice. Keep the jaw surfaces clean and undamaged, because worn jaws hold thread less reliably.

## Supplier

Jain Scientific Equipments Private Limited trades as ScienceLab and supplies the Mixer Artery Forceps from Ambala, Haryana, as part of the vessel-ligation instruments it sources for theatres, training centres, dealers and importers.

Company details are on our [profile page](#).

## Related Surgical Instruments

Instruments

- [Roche](#)
- [vesse](#)

**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana

that clamps



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803