



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Unstirred Thermostatic Water Bath, 22 Litre

Product Code :

SLE/WB/16



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Unstirred Thermostatic Water Bath, 22 Litre (Product Code: SLE/WB/16) is a large analogue laboratory water bath used to warm and incubate many tubes, flasks or bottles at once, from 5°C above ambient to 100°C. Its wide 500 x 300 mm tank opening lets several racks sit side by side, so a whole class or a busy routine lab can share one bath. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, as part of its lab water bath range in India.

Specifications

Parameter	Detail
Capacity	22 litre
Internal dimensions (W x D x H)	500 x 300 x 150 mm
External dimensions (W x D x H)	538 x 332 x 290 mm
Temperature range	Ambient +5°C to 100°C
Sensitivity / stability	±0.5°C at 45°C
Type	Analogue, unstirred, for general use
Operating note	Will operate at 100°C but not designed for continuous boiling; lid recommended above 60°C
Power supply	Mains 230 V AC (frequency listed as 50–56 Hz)
Optional supply	120 V AC, 50–60 Hz models
Electrical safety class	Class 1
Product code	SLE/WB/16
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

The listing for this analogue range gives the following sizes (internal / external dimensions, W x D x H):

- 2.5 litre: 230 x 130 x 150 mm / 270 x 190 x 240 mm
- 4 litre: 300 x 150 x 150 mm / 332 x 185 x 290 mm
- 8 litre: 300 x 240 x 150 mm / 332 x 270 x 290 mm
- 14 litre: 325 x 300 x 150 mm / 361 x 332 x 290 mm
- 22 litre (this model): 500 x 300 x 150 mm / 538 x 332 x 290 mm
- 28 litre: 500 x 300 x 200 mm / 538 x 332 x 290 mm

What Is the 22 Litre Water Bath Used For?

The 22 litre bath is built for volume rather than depth: its 150 mm deep tank has a long 500 mm working length, so rows of test tube racks, culture bottles or flasks can be warmed together.

Typical uses include warming media and buffers for cell culture, pre-warming solutions, and incubating pre-

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803