



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Burette Clamp, Die Cast – Double Clamp

Product Code :

SLE/GBu/11

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Burette Clamp, Die Cast – Double Clamp (Product Code: SLE/GBu/11) is a die-cast retort clamp with a built-in boss head and two spring loaded burette jaws, so a single fitting mounts two burettes side by side on a retort stand rod. The jaws are furnished with rubber grips that close quickly around the glass and hold it steady without scratching, and the metal surfaces carry a high quality finish. Chemistry teachers and analysts looking for a dependable burette clamp for India's school and college titration benches use it for volumetric analysis practicals. It is supplied as a pack of 10 and made by a manufacturer of chemistry lab equipment in India, **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala.

Specifications

Parameter	Detail
Type	Die-cast retort clamp with double burette clamp
Boss head	Built-in
Jaws	Spring loaded, furnished with rubber grips for quick and stable clamping
Finish	High quality finish to metal surfaces
Supplied as	Pack of 10

Typical Use

Slide the built-in boss head onto the rod of a retort stand and tighten it at a height that leaves room for a conical flask on the base. Squeeze a spring loaded jaw open, set the burette in it, and let go: the rubber grips close on the glass and hold it upright. With two jaws on one clamp, students can run an acid burette and an alkali burette side by side for back titrations, or two partners can titrate at the same stand, which halves the number of stands needed on a crowded bench. Whichever retort or burette stand your lab uses, check that each burette hangs vertically before taking the first reading.

Keep the rubber grips clean and free of reagents, as spilled acid or alkali can harden them over time. Wipe the metal after practicals and store clamps dry. The pack of 10 is enough to equip a typical titration class with spares.

Manufacturer

Jain Scientific Equipments Private Limited makes this clamp. The company works under the ScienceLab name from Ambala, Haryana, producing stands, clamps and other apparatus for chemistry and physics teaching. Details are on the [ScienceLab profile page](#).


Related Chemistry Lab Equipment

Items that make up a titration station with this clamp:

- [Burette Clamp, Die Cast – Single Clamp](#) — the single-jaw version for one burette per stand.

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt, Harvna 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803