



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Alcohol for Burning

Product Code :

SCL309/055

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Alcohol for Burning (Product Code: SCL309/055) is a denatured, technical-grade alcohol used as fuel for spirit lamps and alcohol burners in school laboratories. With a strength of 92–96 degrees, it is a standard heat source for practical work in rooms that have no gas supply. Alcohol for Burning is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a lab chemicals supplier in India for schools and colleges.

Specifications

Parameter	Detail
Type	Denatured technical alcohol
Strength	92–96 degrees
Grade	Technical
Intended use	Burning fuel; standard for school laboratories
Shelf life	Shelf life depends on the chemical, its grade and storage conditions; ask us for the expiry date of the batch you are quoted.
Product code	SCL309/055
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is Alcohol for Burning Used For?

Poured into a spirit lamp, this alcohol gives a small, steady flame for warming test tubes, heating water on a tripod and gauze, melting wax in simple heat experiments and flaming an inoculating loop in microbiology. Schools without piped gas rely on it as their everyday bench heat source, and it is equally useful as a backup when gas burners are unavailable.

The liquid is flammable, so fill lamps only when they are cold and extinguished, use a small funnel to avoid drips on the outside of the lamp, and wipe up any spill before relighting. Keep the container tightly closed in a cool, ventilated chemical store away from open flames, and label it clearly: denatured alcohol is not for drinking or internal use.

Supplier

Alcohol for Burning is supplied by Jain Scientific Equipments Private Limited, which trades as ScienceLab from Ambala, Haryana. The company sources laboratory chemicals and consumables so that schools and colleges can order fuel together with the lamps and heating apparatus it goes with. Details about the business are on our [company profile](#).

Related Laboratory Equipment

Apparatus used with burning alcohol at a heating station:

- [Alcohol Lamp](#) — the lamp this fuel is poured into

- [Glass](#)

- [Funne](#)

- [Kamje](#)

- [flame](#)

- [Denat](#)

Website: <https://scilabexport.com>, Email: sales@scilabexport.com over the

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803