



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Filter Paper 1

**Product Code :**

SCL328/064

Image not found or type unknown



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Filter Paper 1** (Product Code: SCL328/064) is a circular laboratory filter paper, 125 mm in diameter, used to separate solids from liquids and to clarify solutions in school, college and applied laboratories. It is widely used in education for teaching simple qualitative analytical separations, and its applications extend to the food industry, agriculture, air pollution monitoring and gas detection. Filter Paper 1 is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a lab consumables supplier in India.

## Specifications

| Parameter         | Detail                                                              |
|-------------------|---------------------------------------------------------------------|
| Form              | Circular filter paper                                               |
| Diameter          | 125 mm                                                              |
| Typical use       | Clarifying liquids; wide range of laboratory applications           |
| Application areas | Food industry, agriculture, air pollution monitoring, gas detection |
| Education         | Teaching simple qualitative analytical separations                  |
| Product code      | SCL328/064                                                          |
| Supplied by       | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India      |

## What Is Filter Paper 1 Used For?

In the classroom, Filter Paper 1 is the sheet that makes filtration visible: chalk or sand stays on the paper while clear water runs through, a coloured precipitate is collected after two solutions are mixed, and in salt analysis the residue and filtrate are tested separately. Beyond teaching, the same paper clarifies liquid samples before they are examined, which is why it appears in food, agricultural and environmental testing work.

For a plain cone, fold the circle in half and then in quarters, open one side and seat it in the funnel, wetting it slightly so it lies flat against the glass. A fluted, fan-style fold gives a larger surface and faster flow when the solid only needs to be removed. Store packs flat and dry, away from chemical fumes that could contaminate the paper.

## Supplier

Filter Paper 1 is sourced and supplied by Jain Scientific Equipments Private Limited, the company that trades as ScienceLab from Ambala, Haryana. Schools, colleges, testing laboratories and dealers can order it together with funnels, flasks and other bench items in a single enquiry, for delivery across India and international markets. Learn more on our [company profile](#).

## Related Laboratory Equipment

Items used with filter paper during filtration:

- [Funnel](#)
- [Conical flask](#)
- [Glass](#)

Website: <https://scilabexport.com>, Email: [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana

trate  
out



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803