



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Havana Porcelain

Product Code :

SCL346/073

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Havana Porcelain (Product Code: SCL346/073) is a porcelain mortar, 80 mm in diameter, accompanied by a 100 mm porcelain press (pestle), used to grind and mix small quantities of solids in school, college and pharmacy laboratories. Its compact size suits crushing crystals, chalk or leaves for a single practical, and it is a familiar item on pharmacy lab equipment lists in India. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

Specifications

Parameter	Detail
Material	Porcelain
Mortar diameter	80 mm
Press (pestle)	Porcelain, 100 mm long, supplied with the mortar
Product code	SCL346/073
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the Porcelain Mortar and Pestle Used For?

Grinding increases the surface area of a solid so that it dissolves, reacts or releases its contents faster. Students use the mortar and press to powder copper sulphate crystals before dissolving them, crush chalk or marble chips for reaction-rate practicals, and grind leaves with a little solvent to extract pigments for paper chromatography.

Grind with a circular pressing motion rather than pounding, and fill the 80 mm bowl only partly so that material does not spill over the rim. Clean the set straight after use; grinding a pinch of clean sand in the bowl helps loosen stubborn residues before washing.

OEM Manufacturer

ScienceLab, the trading name of Jain Scientific Equipments Private Limited, is the OEM for this porcelain mortar and press set, which is made at the company's own licensed factory at Saha, Ambala, Haryana. Buying from the maker means replacement presses or extra mortars for a growing class can be requested from the same source; just ask us. Schools, colleges, dealers and importers are supplied factory-direct, as described on the [ScienceLab company profile](#).

Related Laboratory Equipment

Grinding is usually one step in a longer practical; these items cover the steps before and after it:

- [Porcelain Mortar with Crusher](#) — a comparable porcelain mortar listing.
- [Pestle And Mortar](#) — another mortar option for biology practicals.
- [Spoon Spatula](#) — a porcelain spatula to transfer the ground powder.
- [Porcelain Mortar](#) — a comparable porcelain mortar listing.
- [Plastic Mortar and Pestle](#) — a comparable plastic mortar listing.

Website: <https://scilabexport.com>, Email: sales@scilabexport.com

Find more in [ScienceLab's Laboratory Equipment](#) or [ScienceLab's Laboratory Equipment](#) at 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803