



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Mariot Bottle (for distilled water) 2.5 liter

Product Code :

SCI-603/105

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Mariot Bottle (for distilled water) 2.5 liter (Product Code: SCI-603/105) is a transparent glass storage bottle fitted with a tap, used to hold a working supply of distilled water and run it off in controlled amounts in chemistry, biology and teaching laboratories. Lab assistants rely on it wherever solutions are made up or glassware is rinsed, and the clear glass shows when the 2.5 litre bottle needs refilling. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

Specifications

Parameter	Detail
Capacity	2.5 litre
Dispensing	Fitted with a tap
Material	Glass
Colour	Transparent
Chemical resistance	Chemical resistant
Intended contents	Distilled water
Product code	SCI-603/105
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the 2.5 Litre Mariot Bottle Used For?

Distilled water turns up in nearly every practical: dissolving salts for standard solutions, diluting acids and indicators, making a volumetric flask up to the mark and giving washed glassware a final rinse so tap-water minerals do not affect results. This bottle keeps that water in one clean, clearly identified vessel, and the tap lets a technician fill wash bottles or measuring cylinders without pouring from the neck. As part of school chemistry lab equipment it usually sits in the preparation room or at the end of a bench block.

Because the glass is transparent, cloudiness, floating particles or any growth inside can be spotted before the water is used. Empty and rinse the bottle between fills, run a little water through the tap after cleaning to flush the outlet, and check the tap for drips before leaving the bottle full overnight. Its chemical resistance means it can go through the same cleaning routine as the rest of your laboratory glassware.

What to Specify When Ordering

- Bottle: 2.5 litre with tap, product code SCI-603/105 (a similar distilled-water bottle is listed separately; quote the code you want).
- Quantity: typically one bottle per preparation area or bench block.
- Flow set-up: tell us if you need a constant-flow stopper and tube arrangement as the listing
- Acces

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

OEM Mar 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803