



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Chemical Beaker (Bekera) 100ml

**Product Code :**

SCI-139/042

Image not found or type unknown





Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Chemical Beaker (Bekera) 100ml** (Product Code: SCI-139/042) is a low-form, graduated borosilicate glass beaker used for everyday mixing, dissolving and heating at the student bench in school and college chemistry labs. Its 100 ml capacity makes it the general-purpose size for pair or group work, holding enough solution for a complete reaction while staying light and easy to handle. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

## Specifications

| Parameter           | Detail                                                |
|---------------------|-------------------------------------------------------|
| Capacity            | 100 ml                                                |
| Material            | Borosilicate glass                                    |
| Colour              | Transparent                                           |
| Shape               | Low form                                              |
| Volume markings     | Graduated                                             |
| Chemical resistance | Chemical resistant                                    |
| Heat resistance     | Heat resistant                                        |
| Product code        | SCI-139/042                                           |
| Manufacturer        | Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM |
| Country of origin   | India (made in Ambala, Haryana)                       |

## What Is the 100 ml Chemical Beaker Used For?

Typical jobs for a 100 ml beaker include dissolving a weighed solid in water to make a solution, carrying out a neutralisation with an indicator, collecting the filtrate that drips from a funnel and warming a solution on a tripod and gauze to speed up dissolving. The transparent glass lets students see colour changes, bubbles and precipitates clearly.

Pour from the beaker along a glass rod to avoid drips, keep the liquid well below the rim when stirring, and label beakers that hold different solutions. After use, rinse it with tap water and then distilled water so the next class starts with clean glass. It is one of the core sizes when putting together a beaker set for a school lab in India.

## OEM Manufacturer

Jain Scientific Equipments Private Limited, which does business as ScienceLab, is the original equipment manufacturer (OEM) of the Chemical Beaker (Bekera) 100ml. Production takes place at the company's licensed factory at Saha, Ambala, Haryana, India, so schools and dealers buy from the maker rather than from a trading intermediary. Institutions planning glassware for a new session can ask the manufacturer directly about replacement quantities and repeat orders.

Company based

**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

Related L

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803