



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Dishwashing Brush 200mm

**Product Code :**

SCI-359/146

Image not found or type unknown



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Dishwashing Brush 200mm** (Product Code: SCI-359/146) is a laboratory cleaning brush used to scrub beakers, flasks, dishes and other lab ware after experiments in school, college and research laboratories. It is 200 mm long with a 35 mm brush diameter, a size that reaches into the body of wider vessels, and it is listed as chemical resistant. The brush is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, as one of the everyday lab consumables in its range.

## Specifications

| Parameter           | Detail                                                         |
|---------------------|----------------------------------------------------------------|
| Purpose             | Washing laboratory dishes and glassware                        |
| Brush length        | 200 mm                                                         |
| Diameter            | 35 mm                                                          |
| Chemical resistance | Yes (as listed)                                                |
| Usage               | Laboratory                                                     |
| Product code        | SCI-359/146                                                    |
| Supplied by         | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India |

## What Is a Lab Dishwashing Brush Used For?

Residues left in glassware — dried precipitate, grease, stain or a film of reagent — can spoil the next experiment, so every lab needs a dependable brush at the sink. The 35 mm head of this 200 mm brush suits wide-mouthed beakers, evaporating dishes, Petri dishes and the bodies of conical flasks, where a narrow tube brush would take too long.

Rinse the vessel first, scrub with warm water and a laboratory detergent, then rinse thoroughly with tap water and finish with distilled water before drying on a rack. Use a gentle action on graduated ware so the markings are not worn away, and rinse the brush itself after use so chemical residues are not carried from one vessel to the next.

## What to Specify When Ordering

- Quantity — for example, one per sink or per student group.
- Whether the 200 mm length and 35 mm diameter suit your vessels, or whether narrower tube brushes are also needed.
- Any drying racks or other cleaning items to be quoted together.

## Supplier

Jain Scientific Equipments Private Limited, trading as ScienceLab, supplies this dishwashing brush from Ambala, Haryana, alongside the glassware it is used to clean. Schools and colleges setting up a wash area, as well as dealers and importers stocking lab consumables, can include it in a single order.

The Science Lab is a leading supplier of laboratory equipment and consumables. **Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com) Company background: 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803