



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Burette Clamp - Double

**Product Code :**

SLE/P/26





Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Burette Clamp - Double** (Product Code: SLE/P/26) is a stand-mounted clamp used to hold one or two burettes upright during titrations in school and college chemistry laboratories. It fits stand rods of 10 to 14 mm diameter, and rubber grips keep the burettes from slipping while students read the scale. The clamp is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, the OEM manufacturer, at its own factory in Ambala, India.

## Specifications

| Parameter             | Detail                                                |
|-----------------------|-------------------------------------------------------|
| Burettes held         | One or two                                            |
| Suitable rod diameter | 10 to 14 mm                                           |
| Grip                  | Fitted with rubber grips to prevent slipping          |
| Product code          | SLE/P/26                                              |
| Manufacturer          | Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM |
| Country of origin     | India (made in Ambala, Haryana)                       |

## What Is the Double Burette Clamp Used For?

The clamp attaches to the vertical rod of a retort stand or burette stand and holds burettes vertically above a conical flask. With two positions, a pair of students can share one stand, or one burette can carry the titrant while the second holds another solution, such as a standard for a check titration. Both scales stay at the same height, which makes side-by-side readings easier.

Before fitting, confirm that the stand rod is between 10 and 14 mm thick. Tighten the clamp so each burette hangs truly vertical, then lower it until the tip sits just inside the neck of the flask. The rubber grips hold the glass firmly, so the burette does not slide down as the stopcock is turned. When schools plan chemistry lab equipment for class titrations, a double clamp lets each stand serve two users.

## What to Specify When Ordering

- The rod diameter of your existing stands (this clamp suits 10 to 14 mm rods).
- Whether you need this double clamp, the single version, or a mix of both.
- The number of clamps required per laboratory.
- Whether stands and burettes should be quoted with the clamps.

## OEM Manufacturer

ScienceLab is the trading name of Jain Scientific Equipments Private Limited, the OEM (original equipment manufacturer) of the Burette Clamp - Double. The clamp is made and assembled at the company's own licensed factory at Saha, Ambala, Haryana, and goes factory-direct to schools, colleges, institutions, dealers and importers. If you need replacement rubber grips or extra clamps later, ask us and the manufacturer will handle the request. More about the company is on the [Sci](#)

**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

## Related S

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803