



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Stirrer Bar Retriever

Product Code :

SCI-274/044

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Stirrer Bar Retriever (Product Code: SCI-274/044) is a magnetic pick-up rod used to lift magnetic stirrer bars out of flasks, bottles and beakers in chemistry and biology laboratories. Its strong permanent magnet sits in a PTFE capsule on a screw-in PTFE shaft 8 mm in diameter and 350 mm long, slim enough for narrow-necked vessels, so the bar comes out without tipping out the contents — an everyday companion to any magnetic stirrer. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

Specifications

Parameter	Detail
Magnet	Strong permanent magnet
Magnet capsule	PTFE
Shaft	Screw-in PTFE shaft
Diameter	8 mm
Length	350 mm
Function	Removes stirrer bars from a vessel without tipping out the contents
Especially useful for	Vessels with a narrow opening
Product code	SCI-274/044
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Stirrer Bar Retriever Used For?

Once a solution has been mixed, the stirrer bar must come out before the liquid is made up to the mark, transferred or stored. Pouring risks losing the bar or spilling the sample, and fishing with forceps is awkward in a narrow neck. Lower the retriever until the magnet nears the bar, let the bar jump on to the PTFE tip, and lift it straight out. The 350 mm reach suits volumetric flasks, reagent bottles and round-bottom flasks alike.

PTFE resists most laboratory chemicals and has a non-stick surface, so the retriever can go into acidic, alkaline or solvent-based mixtures and be rinsed clean afterwards. When stirring on a hot plate magnetic stirrer, let the liquid cool before retrieving the bar, and keep the magnet away from electronic balances when it is stored.

Supplier

Jain Scientific Equipments Private Limited, which trades as ScienceLab in Ambala, Haryana, supplies the Stirrer Bar Retriever to school, college, research and quality-control laboratories and to dealers in India and international markets. It is often added to the same enquiry as magnetic stirrers, stirrer bars and glassware. Company information is on the [ScienceLab profile page](#).

Related L

Items used with **Website:** <https://scilabexport.com>, **Email:** sales@scilabexport.com

• **Magn** 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803