



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Petroleum Jelly White

Product Code :

SCL-913/199

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Petroleum Jelly White (Product Code: SCL-913/199) is a soft white semi-solid jelly used to seal, lubricate and protect in school and college laboratories, from bell-jar rims to glass stopcocks. Each 500 gm pack is paraben and phthalate-free, and the product listing also describes it as skin-care jelly for dry skin and minor cuts. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a lab consumables supplier in India for institutional laboratories.

Specifications

Parameter	Detail
Product	Petroleum Jelly White
Capacity	500 gm
Formulation	Paraben and phthalate-free
Skin-care properties (as listed)	Effective in healing dry skin, protecting minor cuts and reducing the appearance of fine lines
Feel on skin	Melts into the skin, making it soft, smooth and visibly healthy
Shelf life	Shelf life depends on the chemical, its grade and storage conditions; ask us for the expiry date of the batch you are quoted.
Product code	SCL-913/199
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is Petroleum Jelly Used For?

A thin smear of petroleum jelly around the ground rim of a bell jar or the flange of a desiccator lid creates an airtight seal, which is essential in experiments on plant respiration, photosynthesis and drying of samples. It is also used to lubricate glass stopcocks and to seal the joints of simple gas apparatus. In biology practicals, students coat the upper or lower surface of leaves with the jelly to block the stomata and compare water loss, and in rusting experiments an iron nail greased with petroleum jelly shows how a barrier layer keeps out air and moisture.

Apply only a thin, even film with a spatula or clean fingertip; excess jelly can creep into apparatus and contaminate samples. Wipe off old grease with a dry cloth before re-applying. Its skin-softening character makes it handy in the lab for chapped hands, although it is not a substitute for gloves when handling chemicals. A pack sits comfortably among the biology lab equipment in India that schools stock for plant physiology work.

Supplier

Petroleum Jelly White is supplied by **Jain Scientific Equipments Private Limited**, trading as ScienceLab from Ambala, Haryana, India, to schools, colleges and institutions in India and international markets. Purchase departments setting the supplier up in their vendor records can note GSTIN 06AAECJ8618A1ZT and Legal Entity Identifier (LEI) 3358004HRD1JIVXKYX49. See the [company](#)

Related C

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803