



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Brush Test Tube

Product Code :

SLE/CH/46



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Brush Test Tube (Product Code: SLE/CH/46) is a hand brush for scrubbing test tubes, made with nylon bristles set on an iron wire stem. It measures 260 mm overall with a bristle head of 90x 45 mm, a compact size that is easy to handle at a student sink. Schools, colleges and routine testing labs use it to remove precipitates, dried solutions and scorch marks after qualitative tests. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, a lab consumables supplier in India that also serves overseas institutions.

Specifications

Parameter	Detail
Type	Test tube cleaning brush
Bristles	Nylon
Stem	Iron wire
Overall length	260 mm
Head size	90 x 45 mm

Typical Use

Salt analysis, flame tests and simple reactions leave test tubes coated with precipitates, oily films or charred residue. Empty the tube into the correct waste container, add a little water and detergent, then twist the brush in and out so the nylon bristles sweep the rounded bottom as well as the walls. Rinse with tap water and finish with distilled water before standing the tube upside down to dry.

Push the brush in gently and do not ram it against the bottom of the tube, because the wire stem can crack thin glass. Stubborn deposits usually come away more easily after a short soak rather than hard scrubbing. For a class set, one brush per bench keeps the clean-up at the end of a practical quick and orderly.

Supplier

ScienceLab is the trading name of Jain Scientific Equipments Private Limited, based in Ambala, Haryana. Alongside the teaching apparatus it manufactures, the company supplies everyday consumables such as cleaning brushes so that institutions can source a complete laboratory order from one place. Details are on the [company profile](#).

Related Chemistry Lab Equipment

Useful companions for washing, handling and drying test tubes:

- [Test Tube](#) — borosilicate test tubes that this brush keeps ready for the next practical.
- [Test Tube Rack With Drying Pegs](#) — lets washed tubes drain on pegs beside the rack holes.
- [Drying Rack](#) — for drying test tubes, flasks, beakers, etc.
- [Test Tube](#) — borosilicate test tubes that this brush keeps ready for the next practical.
- [Brush](#) — for cleaning test tubes, flasks, beakers, etc.

Website: <https://scilabexport.com>, Email: sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803