

Email: sales@scilabexport.com

Phone: +91-7082934803

Product Name:

Product Code:

Portable ECG Recorder w/access

SLE/MDE/13





Email: sales@scilabexport.com

Phone: +91-7082934803

Description:

DESCRIPTION:

Electrocardiogram (ECG) digital monitor and recorder, 12-leads detection, multi-channel recording, portable, AC and battery powered, with printer and accessories

SPECIFICATIONS:

ECG analysis and full interpretation (rhythm and events), real-time and manual ECG diagnostic monitor and recorder, with printer 3- or 6-channels recording, depending on the available model

ECG features:

Common Mode Rejection Ratio (CMRR) > 89 dB frequency response (minimum guaranteed) 0.05 - 150 Hz accuracy of input signal reproduction ± 5 % 12 standard derivations aVR, aVL, aVF, I, II, III and V1 – V6 simultaneous 12-leads acquisition adjustable ECG acquisition and visualisation modes ECG gain 2.5, 5, 10, 20 mm/mV and Auto, accuracy ± 5 % Adjustable ECG sweep (trace speed in mm/s) HR range 30 – 300 bpm with rhythm analysis calibration signal 1 mV ± 5 % input impedance > 2.5 MOhm or > 50 MOhm, depending on the available model pacemaker rhythm detection auto-zero feature

selectable/adjustable filters: baseline correction and Baseline Drift Removal (BDR); AC (mains frequency suppression); Muscle Artefact (EMG) defibrillation and AC fluctuation protection, baseline recovery < 5 s after defibrillation

internal noise level < 15 microV (p-p) current leakage to the patient < 10 microA

Built-in printer:

paper speed (minimum guaranteed) 5, 12.5, 25 and 50 mm/s, accuracy ± 5 % compatible with roll and Z-folded paper (80 mm 3-channel, 110 mm 6-channel) thermal printing head automatic and manual print (selectable modes, ECG traces and reports) 3 or 6 channel recording, depending on the available model

Integrated keyboard

Display, backlit colour TFT LCD, 5 – 8 inch diagonal depending on the available model Visualisation of ECG channels, patient data, date/time, settings, analysis/interpretation, battery status, system messages and alarms

Patient data input

Internal storage for patient data, exams, analysis, interpretations and settings

Data and report export feature

Data/USB and network interface

Adjustable visualisation settings with selectable ECG channels and parameters

and system

Continuous

Audio/visual

Website: https://scilabexport.com, Email: sales@scilabexport.com

er supply

1226/1-5, Bengali Mohalla, Science Market I 133001, Harvana



 ${\bf Email: sales@scilabexport.com}$

Phone: +91-7082934803