



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Spot Check Pulse Oximeter With  
Accessories

**Product Code :**

SLE/MDE/17



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

### DESCRIPTION:

Robust shock resistant design allow use in demanding environments

Suitable for professional clinical use, all patient categories neonate, infant, adult

Ultra-compact pulse oximeter integrated into finger/toe clip

Accommodates finger/toe thicknesses 8 to 25mm directly or through the use of a separate sensor which can be connected to the device.

Spot-check of arterial blood oxygen saturation (SpO<sub>2</sub>) and heart rate(HR)

### Measuring ranges:

**SpO<sub>2</sub> accuracy within:** +/-3%

**Pulse rate detection:** 30 to 240 bpm

**Pulse rate resolution:** 1 bpm)

SpO<sub>2</sub> detection within the range of: 70 to 99%

**Pulse rate accuracy:** within: +/-3%

Display shows SpO<sub>2</sub>, pulse rate, signal strength, sensor error or disconnect and low battery.

SpO<sub>2</sub> resolution: 1% or less

Enclosure protection or better

Suitable for cleaning and disinfection with hospital-grade cleaning products.v

Batteries must allow at least 2500 spot-checks, calculated at 30 seconds per spot-check, or at least 21 hours of operation.

The unit is equipped with an automatic power off feature.

Power requirements two 1.5V AAA batteries

### INSTRUCTIONS FOR USE

Pulse oximeter must be used and maintained properly according manufacturer's instructions.

For spot-check of percentage arterial oxygen saturation (SpO<sub>2</sub>, %) and pulse rate (HR, bpm) of all patient categories. Non-invasive monitor using a clip-on sensor placed on a finger or toe. The spot-check pulse oximeter is ultra-compact and fits easily pocket or carry-on bag.

### ENVIRONMENTAL CONDITIONS

**Storage conditions temperature:** -20 °C to 60°C / RH: 10% - 95%

**Operating conditions temperature:** 5 to 40°C / RH: 30 to 80%

Contact Science Lab for your Educational School Science Lab Equipments. We are best technical educational equipments manufacturers, technical educational equipments manufacturer, technical educational equipments supplier, technical educational instruments, technical educational instruments exporter, technical educational instruments manufacturer, vocational training lab instruments exporter