



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Rubber Stoppers without Holes Of Different
Diameters

Product Code :

SCI-204/024

Image not found or type unknown





Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Rubber Stoppers without Holes Of Different Diameters (Product Code: SCI-204/024) are solid, round rubber stoppers used to close test tubes, small flasks and bottles in school, college and research laboratories. Having no bore, they simply shut the opening, and they come in bottom diameters of 11 mm, 13 mm, 14 mm and 15 mm for storing, shaking and protecting samples. They are supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a lab consumables supplier in India serving schools and institutions.

Specifications

Parameter	Detail
Type	Solid stopper, without hole
Material	Rubber
Shape	Round
Chemical resistance	Chemical resistant
Bottom diameters available	11 mm, 13 mm, 14 mm, 15 mm
Use	Laboratory
Product code	SCI-204/024
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Are Solid Rubber Stoppers Used For?

Solid stoppers do the quiet jobs of a practical class. Students stopper a test tube before shaking it to dissolve a solid or compare solubility, close a tube of collected gas until it is ready for the splint test, and seal small reagent or sample bottles so that dust stays out and evaporation is reduced. They also blank off a spare neck while apparatus is being assembled, which is why technicians keep a box of each size in the preparation room.

Never close a tube or flask with a solid stopper while it is being heated or while a reaction is still giving off gas: pressure builds up and can push the stopper out or crack the glass. Press the stopper in with a slight twist rather than a hard push, and ease it out the same way. If a set-up later needs a tube through the stopper, a hole can be bored with a wetted cork borer.

What to Specify When Ordering

- Bottom diameter(s): 11 mm, 13 mm, 14 mm or 15 mm.
- Quantity of each size for every lab or bench.
- The glassware to be closed, so sizes can be checked against its neck.
- Whether the one-hole version is also needed for delivery tubes and thermometers.

Supplier

Jain Scientific Equipments Private Limited, trading as ScienceLab, supplies these solid rubber stoppers from its Ambala factory to schools, colleges and research laboratories in India and internationally. **Website:** <https://scilabexport.com>, **Email:** sales@scilabexport.com Often order the solid and liquid stoppers, test tubes, flasks, beakers, etc. Registration, 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803