



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Measuring Cylinder 1000 ml

Product Code :

SCI-325/053

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Measuring Cylinder 1000 ml (Product Code: SCI-325/053) is a one-litre graduated borosilicate glass cylinder used to measure bulk liquid volumes and to run tall-column experiments in school, college and research laboratories. With 10 ml divisions, it lets teachers show exactly what a litre looks like and gives physics and geography classes a tall column for viscosity and sedimentation work. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India, for school lab setups at home and abroad.

Specifications

Parameter	Detail
Capacity	1000 ml (one litre)
Graduation interval	10 ml (listed as "glass separation rate")
Material	Borosilicate glass
Heat resistance	Heat resistant
Chemical resistance	Chemical resistant
Marking	Graduated
Product code	SCI-325/053
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is a 1000 ml Measuring Cylinder Used For?

Filled to the top mark, it shows students that one litre equals one cubic decimetre, a simple way to connect units with a real volume. Prep rooms use it to measure distilled water for large batches of solution. In physics it becomes the column for the falling-ball viscosity experiment, and in geography or agriculture lessons a shaken soil-and-water mix settles in it into visible sand, silt and clay layers.

When full it holds about a kilogram of water on top of the weight of the glass, so lift it with both hands, one under the base, or better still fill and empty it where it stands. Clean glycerine or oil out with warm soapy water, and keep the cylinder back from the bench edge when not in use.

What to Specify When Ordering

- Number of litre cylinders for demonstration benches and the prep room
- Smaller capacities from the series you want in the same quotation
- Whether plastic cylinders or cleaning brushes should be quoted alongside

OEM Manufacturer

Jain Scientific Equipments Private Limited trades under the ScienceLab name and manufactures this litre cylinder. We are the maker and sell direct, source, and glassware —

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803