



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803


**Product Name :**

Marked (tarred) Pots 1000 ml

**Product Code :**

SCI-638/090

Image not found or type unknown



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



**Marked (tarred) Pots 1000 ml** (Product Code: SCI-638/090) is a litre-size transparent borosilicate glass tare vessel used to make up and hold bulk quantities of solutions in school, college and institutional laboratories. It closes with a sanded cap and is chosen when a full day of practicals draws on one shared stock. This borosilicate glassware is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

| Parameter         | Detail                                                |
|-------------------|-------------------------------------------------------|
| Capacity          | 1000 ml                                               |
| Material          | Borosilicate glass                                    |
| Colour            | Transparent                                           |
| Design            | Tare vessel with sanded cap                           |
| Resistance        | Chemical resistant, heat resistant                    |
| Product code      | SCI-638/090                                           |
| Manufacturer      | Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM |
| Country of origin | India (made in Ambala, Haryana)                       |

Laboratories use this pot for litre-scale batches: a technician tares the empty vessel, weighs in the solute and makes the solution up before the first class, so every group draws from the same stock. It also works as a clean reservoir for distilled or rinse water on a busy bench. Because the listing for this size does not mention graduations, measure volumes with a measuring cylinder or volumetric flask.

Borosilicate glass resists heat and most routine chemicals, which is why it is preferred for chemistry lab glassware that is washed and reused for years. Keep the sanded cap and the rim clean and dry, because grit on a ground surface can make a cap stick. Lift the full vessel by its body, never by the cap.

- Capacity: 1000 ml under code SCI-638/090; smaller marked pots are separate listings.
- Graduations: if you need volume marks at this size, say so in your enquiry so we can confirm what applies.
- Quantity: the number of litre-size stock vessels each lab needs.

Jain Scientific Equipments Private Limited operates under the ScienceLab name and is the OEM (original equipment manufacturer) of the Marked (tarred) Pots 1000 ml. The vessel is made at the company's li [REDACTED] from a trader, and ships di [REDACTED] and matching [REDACTED].  
**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803