



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Washing Bottles (Trays)

Product Code :

SCI-303/120

Image not found or type unknown



Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Washing Bottles (Trays) (Product Code: SCI-303/120) is a 500 ml plastic squeeze-and-spray wash bottle used to rinse glassware, wash down precipitates and dispense distilled water or solvents in a controlled stream in chemistry and biology laboratories. Its sealable screw-top helps stop solvents from evaporating, and a 90° nozzle bend lets the user aim the spray at different angles. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer of lab plastic ware, at its own factory in Ambala, India.

Specifications

Parameter	Detail
Material	Plastic
Capacity	500 ml
Type	Squeeze-and-spray bottle
Closure	Sealable screw-top that prevents solvents from evaporating
Nozzle	90° bend for multiple spray angles
Properties	Chemical resistant
Product code	SCI-303/120
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is a 500 ml Wash Bottle Used For?

A wash bottle is the tool for any job that needs a small, directed stream of liquid: the final distilled-water rinse of a beaker, washing a precipitate from the walls of a beaker onto filter paper, transferring every trace of a sample into a volumetric flask, or topping up a solution close to a graduation mark. The bent nozzle reaches into flasks and funnels without tipping the bottle over.

Keep one bottle for each liquid and label it clearly, for example distilled water or ethanol, so contents are never mixed. Squeeze gently for a fine jet, keep the screw-top sealed when the bottle is not in use, and refill from a clean stock container rather than topping up an old batch.

What to Specify When Ordering

- The capacity: this listing is the 500 ml bottle; other plastic washing bottles are separate items.
- How many bottles you need per bench, allowing one for each liquid in use.
- Any funnels or distilled-water containers to include in the quotation.

OEM Manufacturer

Washing Bottles (Trays) are produced by Jain Scientific Equipments Private Limited, trading as ScienceLab, at the company's own factory in Ambala, Haryana. You can order straight from the factory to get the best price. The factory is located at 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana. It can be arranged with

Website: <https://scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803