



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Centrifuge Tubes

**Product Code :**

SCI-222/034

Image not found or type unknown



**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Centrifuge Tubes** (Product Code: SCI-222/034) are sterile 15 ml polypropylene tubes with seal caps used to hold, spin and store liquid samples in biology, microbiology and chemistry laboratories. The optically clear, radiation-sterilised tubes carry silk-screened exterior volume graduations and a large white frosted area for labelling, and a plug built into each cap makes the closure leak-proof. They are supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

## Specifications

| Parameter        | Detail                                                         |
|------------------|----------------------------------------------------------------|
| Capacity         | 15 ml                                                          |
| Material         | Optically clear polypropylene                                  |
| Closure          | Seal cap with built-in plug for a leak-proof fit               |
| Graduations      | Silk-screen exterior volume graduations                        |
| Marking area     | Large white frosted portion for easy writing                   |
| Sterility        | Sterile; radiation sterilised                                  |
| Autoclaving      | Can be autoclaved                                              |
| Composition note | Free from natural rubber and heavy metal                       |
| Product code     | SCI-222/034                                                    |
| Supplied by      | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India |

## What Are 15 ml Centrifuge Tubes Used For?

Tubes of this size are the everyday workhorse of a laboratory centrifuge: separating cells or precipitates from a liquid, pelleting bacteria from broth cultures, clarifying extracts and holding samples between steps. The clear polypropylene lets students see the pellet and the liquid level, while the printed graduations give a quick volume reading without a separate measuring cylinder.

Always balance tubes of equal volume opposite each other in the rotor, check that your centrifuge accepts 15 ml tubes, and write sample details on the frosted panel before filling. The sealed cap helps when tubes are shaken, inverted or stored, which makes them a useful stock item for anyone running a laboratory centrifuge in India or abroad.

## What to Specify When Ordering

- The quantity of 15 ml tubes required.
- The centrifuge model or rotor you use, so compatibility with 15 ml tubes can be confirmed.
- Whether you also need tube racks or pipettes for filling.

## Supplier

Jain Scientific  
centrifuge tu  
vendor regis

**Website:** <https://scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana

the sterile  
s. For  
ed (CIN



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803