



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Electromagnet with Transformer for
Sonometer

Product Code :

SCL-292/044

Image not found or type unknown





Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Electromagnet with Transformer for Sonometer (Product Code: SCL-292/044) is a sonometer accessory used to set a steel sonometer wire vibrating from the AC mains, so students can determine the frequency of the supply in the physics lab. Senior schools and colleges use it for the AC sonometer experiment, and it pairs the electromagnet with its own transformer in one item. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India, for sonometer work in Indian and overseas labs.

Key Details

Parameter	Detail
Item	Electromagnet with transformer
Components	Electromagnet; transformer
For use with	Sonometer
Product code	SCL-292/044
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the Sonometer Electromagnet Used For?

In the AC sonometer experiment the electromagnet sits close to the steel wire between the two bridges, and the transformer feeds its coil from the AC mains. The alternating field pulls on the magnetic wire twice in every cycle, so the wire is driven at twice the mains frequency. Students move a bridge or change the load until a small paper rider on the wire is thrown off, which marks resonance, then calculate the wire's frequency from its vibrating length, tension and mass per unit length and halve it to get the supply frequency.

The method only works with a magnetic wire such as steel; brass or copper wire will not respond to the electromagnet. Keep the pole close to the wire without touching it, switch off at the mains before adjusting the set-up, and avoid leaving the coil energised between readings so it does not warm up unnecessarily.

What to Specify When Ordering

- Your mains supply voltage and frequency, so we can confirm the transformer suits it
- The sonometer it will be used with, for example a simple or 3-wire pattern
- Whether a sonometer, steel wire and slotted weights should be quoted with it

OEM Manufacturer

ScienceLab, the trading name of **Jain Scientific Equipments Private Limited**, is the OEM of this Electromagnet with Transformer for Sonometer and makes it at its own licensed factory at Saha, Ambala, Haryana, India. Buyers can put questions on compatibility with their sonometer, replacement parts and request orders straight to the manufacturer's sales. The company holds MSME registration and exports to over 100 countries. It is an ISO 9001:2015 certified company and ships export orders worldwide. For more information, visit the [ScienceLab website](https://scilabexport.com) or contact us on WhatsApp at +91 7082934803. **Website:** <https://scilabexport.com>, **Email:** sales@scilabexport.com and is on the [ScienceLab website](https://scilabexport.com). **Address:** 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803