Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Drying Rack Glass & Plastic Ware

Product Code :

SLE/CL/54



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Drying Rack Glass & Plastic Ware (Product Code: SLE/CL/54) is an upright, plastic-coated wire rack on which freshly washed beakers, flasks, cylinders, bottles and plastic containers are inverted to drain and air-dry beside the sink. It suits wash-up areas in clinical, school and college laboratories, where it is often ordered together with new glassware from a single lab glassware supplier in India. The rack is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala.

Specifications

Parameter	Detail
Material	Plastic-coated wire
Mounting	Free-standing or wall mount
Dimensions (approx.)	36 x 16 x 65 cm (w x d x h)
Intended for	Drying washed glass and plastic ware

Typical Use

After washing and a final rinse, turn each vessel mouth-down on the rack so that water runs out rather than collecting inside, and hang narrow-necked items such as volumetric flasks and measuring cylinders so they cannot tip over. The plastic coating keeps bare metal away from the glass, which reduces knocks and scratches while wet items are being handled, and it is equally suitable for polypropylene beakers, bottles and sample containers.

Stand the rack on a draining board next to the sink, or fix it to the wall above the sink when bench space is short; water then drips straight into the basin. Keep heavier glassware on the lower part of the rack for stability. Schools and diagnostic labs that buy from a lab plastic ware supplier in India often pair this rack with a wash bottle of distilled water so that the final rinse leaves no mineral spots as the items dry.

Manufacturer

ScienceLab is the brand under which Jain Scientific Equipments Private Limited, based in Ambala, Haryana, makes wash-room, bench and teaching equipment for laboratories. Visit our company profile for background on the business and its product range.

Related Clinical Laboratory Equipment

Items that belong to the same wash-and-dry routine:

- Brush Test Tube — scrubs the inside of tubes and narrow vessels before they go on the rack.
- Wash Bottle LDPE with Narrow Mouth — directs a final rinse of distilled water over cleaned labware.
- Beaker Plastic — everyday plastic beakers that are dried on this rack between uses.
- [illegible]e
- [illegible]