Email : sales@scilabexport.com
Phone: +91-7082934803

**Product Name :**

Vacuum Extractor Bird Anterior Manual Complete Set

**Product Code :**

SLE/HE/02



Email : sales@scilabexport.com
Phone: +91-7082934803

## Description :

**Vacuum Extractor Bird Anterior Manual Complete Set** (Product Code: SLE/HE/02) is a hand-operated vacuum extraction set with Bird-type anterior cups, used to assist vaginal deliveries in hospital delivery rooms and maternity units. It pairs stainless steel anterior cups of 50 and 60 mm with two soft cups and a 0 to -100 kPa (-800 mmHg) vacuum gauge, and all of its parts can be disassembled and autoclaved at 121 °C. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

### Specifications

| Parameter | Detail |
|---|---|
| Application | Vacuum extraction to assist vaginal deliveries |
| Vacuum source | Hand-operated suction pump |
| Vacuum gauge range | 0 to -100 kPa / -800 mmHg |
| Cup type | Bird-type anterior cups, stainless steel, plus soft cups |
| Reprocessing | All parts can be disassembled and are autoclavable at121 °C |
| Maintenance | No specific maintenance required |
| Storage and transport | -20 °C to +50 °C, 20% to 95% humidity |
| Weight and volume (as listed) | 29 kg; 0.73 m³ — confirm packed shipping details on the quotation |
| Product code | SLE/HE/02 |
| Supplied by | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India |

Complete set contents:

- 1 × vacuum gauge (0 to -100 kPa / -800 mmHg) with straight connector, chrome plated
- 1 × vacuum bottle, plastic, 500 ml
- 1 × cover plug for the collection bottle, with 3 holes
- 2 × Bird-type anterior suction cups, stainless steel, 50 mm and 60 mm diameter, rounded edges
- 2 × soft vacuum extractor cups, 50 mm and 60 mm diameter
- 2 × extraction handles, stainless steel
- 1 × angled connector, chrome plated
- 1 × angled connector with screw valve, chrome-plated screw, cap and gasket
- 1 × epoxy-coated basket for the bottle and suction pump
- 2 × transparent silicone tubes, inner/outer diameter 6/11 mm, lengths 50 cm and 150 cm

### What Is a Bird Anterior Vacuum Extractor Used For?

The set is assembled as a closed vacuum line: the chosen cup connects by silicone tube and angled connector to the 3-hole cover plug of the 500 ml collection bottle, the gauge connects throu xy-coate s suit d