Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Hospital Bed Screen On Castors

Product Code :

SLE/HF/02



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Hospital Bed Screen On Castors (Product Code: SLE/HF/02) is a mobile, three-section privacy screen fitted with non-transparent plastic curtains and mounted on swivel castors, so ward staff can wheel it to any bed and shield the patient during examination, dressing changes or personal care. It is basic equipment for general care inpatient wards in hospitals and health centres, and hospital and nursing-college buyers handling laboratory equipment tenders in India often list it alongside beds and bedside cabinets. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)** in Ambala, India.

Specifications

Parameter	Detail
Configuration	3 sections with curtains
Panel size	3 panels of approx. 600 × 2100 mm (w × h) each
Mobility	Swivel castors, approx. 50 mm diameter, positioned in a broad stance for good stability
Frame	Epoxy or enamel coated, or chrome plated, tubular steel
Curtains	Plastic, non-transparent, washable and flame retardant
Material resistance	Resistant to corrosion (tropical environment) and to decontamination with chlorine solution
Construction	Knockdown

Supplied with:

- 1 set of plastic curtains
- Clear instructions and diagrams for assembly
- List of accessories and parts, as well as specific tools for assembly

Typical Use

On an open ward, nurses roll the screen to the bedside and arrange its three panels around the patient to create a private space for physical examination, catheter or wound care, bathing and bedpan use. Because the curtains are opaque and the panels stand about as tall as a doorway, the patient is shielded from the rest of the ward without the need for fixed ceiling tracks, which makes the screen useful in older buildings, temporary wards and outreach clinics. Nursing schools also use it to teach dignity-in-care routines in their practice wards.

The knockdown frame ships flat and is assembled on site using the diagrams provided. Wipe the curtains and frame with a chlorine solution between patients, and check the castors periodically so the screen stays stable when it is moved or nudged.

Manufacturer

The s
from
More