Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Instrument Table SS 2 Trays On Castors

Product Code :

SLE/HF/18



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Instrument Table SS 2 Trays On Castors (Product Code: SLE/HF/18) is a mobile stainlesssteel table with two flat, fixed shelves, used to lay out instruments, dressings and supplies withinreach of the clinician. Operating theatres, minor-procedure rooms, clinics and nursing or dentaltraining labs use it as a wheeled working surface that is quick to wipe down between cases. It ismade by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, one of the lab instruments manufacturers in India that also builds furniture for hospitals and clinical training rooms.

Specifications

Parameter	Detail
Type	Movable instrument table with two flat, non-removable shelves (trays)
Top shelf	No raised edges — an entirely flat work surface
Lower shelf	Fully flat, fixed to the table frame just above the swivel castors
Shelf finish	Smooth finish for easy cleaning; 1.5 mm thickness
Frame and trays	Austenitic stainless steel 18/10; frame 3 cm (outside, across), 1.5 mm thickness
Bracing	Transfer bars connect the lower distal portions of all 4 legs
Castors	Heavy carriage on 4 swivel castors, all 4 with anti-static wheels, 2 with total-lock brake (wheel and rotational lock)
Castor construction	Single mold-on wheel, non-hooded for easy maintenance; sealed bearings in swivel and wheel; ball-bearing swivel; frame/bracket steel or nylon
Castor wheel size	2.5 x 10 cm (w x diameter)
Overall size	90 x 50 x 92 cm (l x w x h)
Carrying capacity	100–120 kg
Environment	High resistance to corrosion (tropical environment)

Typical Use

Before a procedure, a drape is spread over the flat top shelf and the instrument set is opened onto it in the order it will be used; the lower shelf carries spare packs, gauze and solutions so the top stays uncluttered. Because the top has no rim, large trays and drapes lie flat to the edge, and the table can be wheeled close to the operating or examination couch and braked in position.

Beyond theatres, the same table works as a mobile workstation in teaching and research settings — lab equipment for research in India often includes a stainless steel table like this for dissection instruments, specimen trays or demonstration kits. After use, clear both shelves, wash them with a neutral detergent, disinfect, and dry; steel wool and chlorine-based cleaners left standing can mark stainless steel over time.

Man

Jain S is a