Science Lab

Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Non-Pneumatic Anti-Shock Garment

Product Code :

SLE/DDS/014



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Non-Pneumatic Anti-Shock Garment (Product Code: SLE/DDS/014) is a reusable neoprene garment with Velcro fastenings used as a first-aid measure for women with uncontrollable postpartum haemorrhage (PPH) in maternity units, health centres and referral transport. Listed in a small size, it is intended to help control the bleeding, reverse the shock and stabilise the patient for safe transport to comprehensive obstetric care, and it has no metal parts. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

Key Details

Parameter	Detail
Type	Non-pneumatic anti-shock garment (NASG)
Size	Small
Intended use	Women with uncontrollable postpartum haemorrhage (PPH): to control the bleeding, reverse the shock and stabilise the patient for safe transport to a comprehensive obstetric care facility
Material	Lightweight neoprene garment with Velcro
Metal parts	None; described as safe for X-rays and MRIs
Wear	Lightweight, flexible and comfortable for the wearer
Use	Reusable
Product code	SLE/DDS/014
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Anti-Shock Garment Used For?

Postpartum haemorrhage can become critical faster than a woman can reach a facility able to treat its cause. The non-pneumatic anti-shock garment is intended for that gap: it is applied to a woman with uncontrollable PPH to help control the bleeding and reverse the shock, so that she can be stabilised and transported safely to comprehensive obstetric care.

Being non-pneumatic, the garment needs no pump or pressure gauge; the neoprene is closed with Velcro. It has no metal parts, so it is described as safe for X-rays and MRIs, and it is reusable. It should be applied and removed only by health workers trained in its use, following the protocol of their facility or maternal-health programme.

What to Specify When Ordering

- Number of garments required
- Size: this listing is the small size – ask us about any other sizes before ordering
- Other delivery-room items to be quoted together
- Delivery destination

Sup

Jain S ia, suppl nal-