Science Lab

Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Micro Centrifuge Tube, 1.5ml

Product Code :

SLE/P/41



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Micro Centrifuge Tube, 1.5ml (Product Code: SLE/P/41) is a small polypropylene tube with aflat cap used to hold, spin and store microlitre-scale samples in biology, biotechnology andclinical laboratories. Packed 500 to a pack, the tubes are made from robust virgin polypropylenewith clear moulded graduations and a frosted writing area. These lab consumables are suppliedby **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

Specifications

Parameter	Detail
Capacity	1.5 ml
Material	Robust virgin polypropylene
Clarity	Clear, for better visual inspection of the sample
Cap	Flat cap with a bevelled edge for comfortable opening and closing
Cap top	Flat, with a writable surface
Labelling area	Frosted writing area on the tube
Graduations	Clear moulded graduations on the tube
Sealing	Precise lid sealing to keep evaporation low during storage
Pack size	500 tubes
Product code	SLE/P/41
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is a 1.5 ml Micro Centrifuge Tube Used For?

Micro centrifuge tubes are the standard vessel for small-volume work such as DNA and protein preparations, enzyme reactions, spinning down cells or precipitates, and keeping aliquots for later use. Students pipette the sample in, close the cap, label the frosted area or the cap top and then centrifuge, mix or incubate the tube.

Balance the rotor by loading tubes of equal volume opposite each other, and check that your rotor takes 1.5 ml tubes before a run. Most labs treat these tubes as disposable lab equipment and order full packs for a term of practicals, which is why they come in packs of 500.

What to Specify When Ordering

- The number of packs — each pack contains 500 tubes of 1.5 ml.
- Whether sterile tubes are needed; the listing does not state sterility, so ask us.
- The centrifuge or rotor the tubes must fit, so compatibility can be checked.
- Matching items such as a micropipette, tips or a tube rack.

Supplier

Jain S l
tubes
deale m