Email : sales@scilabexport.com
Phone: +91-7082934803

**Product Name :**

Spatula, 140mm

**Product Code :**

SLE/CH/222



Email : sales@scilabexport.com
Phone: +91-7082934803

## Description :

**Spatula, 140mm** (Product Code: SLE/CH/222) is a 140mm stainless steel laboratory spatula with a round bent end at one side and a scoop end at the other. The scoop end carries a small portion of powder or crystals from a reagent bottle, while the bent end is handy for scraping, stirring or breaking up caked solids. Its silver stainless steel finish is easy to wipe clean, making it an everyday tool at the chemical balance or digital weighing balance in school, college and industrial labs. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

### Specifications

| Parameter | Detail |
|---|---|
| Material | Stainless steel |
| Colour | Silver |
| Length | 140mm |
| End 1 | Round bent end |
| End 2 | Scoop end |

### Typical Use

A typical task is weighing out a salt to make up a standard solution: the student tares a weighing boat, uses the scoop end to add solid a little at a time, and taps the spatula gently to add the last few crystals until the target mass is reached. The round bent end is useful for loosening a reagent that has caked inside its bottle, for stirring a small quantity of mixture in a beaker, or for scraping residue from the sides of an evaporating basin.

To avoid cross-contamination, wipe and dry the spatula each time you change chemicals, or keep a dedicated spatula for each frequently used reagent. Never return unused solid to the stock bottle from the spatula; place it in a labelled container instead. Wash with detergent and water after the session and store dry.

### Manufacturer

This spatula is made by Jain Scientific Equipments Private Limited, the company behind the ScienceLab name in Ambala, Haryana. The company is incorporated under CIN U29309HR2020PTC087597 and manufactures at its licensed factory at 449, HSIIDC Industrial Area, Saha, Ambala (Factory Licence AMB-ONLINE-CHD-J-414). Company details are on the profile page.

### Related Chemistry Lab Equipment

Items used with this spatula when weighing and transferring solids:

- Plastic Weighing Boats — receive the solid on the balance pan.
- Balance Double Beam — for weighing the portions the spatula delivers.
- y.
- 
- 