
Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Lens Tissue

Product Code :

SLE/M/121



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Lens Tissue (Product Code: SLE/M/121) is a soft, disposable cleaning paper used to wipe camera, spectrometer and microscope lenses in laboratories and classrooms. Each sheet measures 100 x 150 mm and the tissue comes in a pack of 100, a practical quantity for a microscope lab's routine cleaning. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, a laboratory instruments supplier in India, as a consumable for optics and biology labs.

Specifications

Parameter	Detail
Sheet size	100 x 150 mm
Pack size	Pack of 100
Used for	Cleaning camera, spectrometer and microscope lenses
Product code	SLE/M/121
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is Lens Tissue Used For?

Lens tissue is the everyday wipe for optical glass: microscope eyepieces and objectives inbiology labs, camera lenses on imaging setups, and the lenses of a spectrometer in physics practicals. Ordinary paper tissue and clothing can leave fibres or drag grit across coated surfaces, which is why labs keep proper lens tissue beside every microscope bench.

Use each sheet once. Blow off loose dust first, fold a sheet so your fingers stay away from the wiping area, add a drop of lens cleaning fluid if needed, and wipe gently in a spiral from the centre outward. To lift immersion oil from an objective, draw a fresh strip of tissue across the front lens, then finish with a dry sheet.

Supplier

Jain Scientific Equipments Private Limited, which trades as ScienceLab from Ambala, Haryana, supplies lens tissue as a consumable alongside its microscopes, cameras and optics equipment. Schools, colleges, research labs, dealers and importers in India and international markets can add packs to larger lab orders. More is on the ScienceLab company profile.

Related Microscope Accessories

Other items used with lens tissue:

- Microscope Cleaning Kit — brush, air blower, lens fluid and swabs for a full clean
- Immersion Oil, 25ml — wiped off oil objectives with lens tissue after use
- Microfiber Lens Cloth — reusable cloth for eyepieces and instrument bodies
- Spectrometer — physics instrument whose lenses need the same care

Find

Frec