Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Hazardous Storage Cabinet Large

Product Code :

SLE/HF/48



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Hazardous Storage Cabinet Large (Product Code: SLE/HF/48) is a robust, horizontal two-door steel cabinet for keeping flammable liquids, chemicals, paints, solvents and other dangerous substances locked away from general storage. Made from 0.9mm thick steel and finished inbright yellow with warning signs, it is designed for good structural stability and has a leak proofsump that contains spillages inside the cabinet. Hospitals, laboratories, maintenance departmentsand schools use it as a lab storage cabinet in India wherever hazardous stock must besegregated. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala.

Specifications

Parameter	Detail
Type	Horizontal
Material	Steel, 0.9mm thick
No. of doors	2
No. of shelves	1 (adjustable, perforated, galvanised)
Sump	Leak proof, removable, galvanised
Finish	Bright yellow with warning signs
Dimensions (H x W x D)	1000 x 915 x 457mm
Suitable for	Flammable liquids, chemicals, paints, solvents and dangerous substances

Typical Use

The horizontal, waist-height format (1000mm high and 915mm wide) makes this the model to choose when the cabinet must sit under a window, beside a workbench or in a room where a tall unit would be unstable or block sightlines. Typical contents include solvent bottles in a pathology or histology lab, spirit and disinfectant stocks in a hospital store, paints and thinners in a maintenance workshop, and flammable reagents in a school chemistry prep room. The bright yellow body and warning signs tell staff at a glance that the contents need care.

Use the adjustable perforated shelf to split tall containers from small bottles, and stand leaking or opened containers so that any drips fall into the sump, which can be lifted out for cleaning. As with any chemistry reagent cabinet in India, keep incompatible groups apart — acids and oxidising agents should not share space with flammable solvents — and place the cabinet away from heat sources and exits.

Manufacturer

Hazardous storage cabinets in this range come from **Jain Scientific Equipments Private Limited**, which markets its products under the **ScienceLab** name from Ambala, Haryana, and serves hospitals, laboratories and educational buyers. Visit the company profile to learn more about

Rela