Email : sales@scilabexport.com
Phone: +91-7082934803

Product Name :

Burette Brush

Product Code :

SLE/GBu/10



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Burette Brush (Product Code: SLE/GBu/10) is a long cleaning brush for the narrow barrel of a burette, measuring 24" overall with an 8.5" section of white nylon bristles and a 0.50" brush diameter. Its double thick gauge, wire-wound handle is heavy-duty yet flexible, so it can be pushed down a full-length burette tube. Laboratory technicians and students use it to clear reagent films and residues after titrations. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, a lab consumables supplier in India and international markets.

Specifications

Parameter	Detail
Overall length	24"
Bristle section	8.5", nylon bristles
Brush diameter	0.50"
Bristle colour	White
Handle	Double thick gauge wire-wound, heavy-duty flexible metal
Purpose	Cleaning burettes and similar long, narrow glassware

Typical Use

After a titration session, rinse the burette with tap water, wet the brush in a mild laboratory detergent solution and work it down the barrel from the top with a gentle twisting motion, then draw it back. The long handle reaches the lower graduations without the hand entering the tube, and the flexible wire lets the brush follow the tube rather than jamming. Finish with tap water and a final rinse of distilled water, and let the burette drain inverted.

Do not force the brush into the narrow jet below the stopcock; clean that part by flushing instead. Keep the twisted wire tip away from hard contact with the glass, rinse the brush after each use and hang it to dry. When bristles begin to splay or wear, replace the brush — a routine item for any chemistry consumables supplier in India to stock alongside other glassware cleaning tools.

Supplier

This brush is supplied by **Jain Scientific Equipments Private Limited**, which trades as **ScienceLab** from Ambala, Haryana, India, and supports laboratories with glassware, apparatus and cleaning accessories. For vendor registration forms, the company's GSTIN is 06AAECJ8618A1ZT and its LEI is 3358004HRD1JIVXKYX49. Read more on our company profile.

Related Laboratory Equipment

Other cleaning tools and the glassware this brush is used on:

- Brush Burette — another burette brush option in our chemistry range
- Brush Test Tube — a shorter brush for test tubes